



Cleaning



Tablets



Pumps



Ancillaries



Ice machine sanitising cleaner concentrate



CONCENTRATED

NON HAZARDOUS



10:1

Makes 2.5 litres



Non-corrosive, safe on all metals, aluminium and copper. Safe on all surfaces.

Ice-clean is a concentrated disinfectant, sanitiser and deodoriser, specifically designed for ice machines in bars, restaurants, and food retail. Ice-clean reduces the danger of transmission of many infectious organisms and cross contamination from environmental surfaces. The solution improves system efficacy and keeps ice hygienic. The solution is odourless and food safe.

Ice-clean sanitises ice machines to the highest standards, breaking down the biofilms that allow bacteria to survive. Biofilms that develop in ice machines can contaminate the ice, and reduce efficiency. Cleaning with Ice-clean will deliver a long lasting antibacterial coating to prevent and control biofilms.


Effective to BS EN1276 & EN1650 standards.
NSF Registered for food preparation areas.




Nonfood Compounds D2



Add to reservoir at a dilution of 10:1.



Follow manufacturers guidelines for cleaning cycle.



Rinse with clean water before making ice.

Effective against bacteria, mould, fungi & viruses, including Human Coronavirus, Hepatitis B, Hepatitis A, E-Coli, Salmonella, Herpes, Avian flu, H3N2 flu (Hong Kong), Listeria, MRSA, HIV-1 & HIV-2.

Long lasting antibacterial protection remaining effective long after cleaning.



Name	Ice-clean concentrate
Product code	ICC250
Weight	265 grams
Colour	Clear
Hazards	None
Soils	Dust, Airbourne pollutants, bacteria, mould, slime & biofilms

Ice-clean concentrate is a non-hazardous concentrate requiring no COSHH / REACH controls for transport as storage.

INGREDIENTS;

Aqua, Urea, Benzalkonium chloride, Active ingredients <1% by volume.

CAUTION; Keep out of reach of children. Do not swallow. Avoid eye contact. In case of contact with eyes, rinse with water, if swallowed, wash mouth out with water and seek medical advice if symptoms arise or persist.

STORAGE; Store upright at room temperature. Do not allow to freeze.

0844 330 3643 info@gel-clear.co.uk gel-clear.co.uk